

## Recommended Steps:


# Putting On Personal Protective Equipment (PPE)

### 1. Perform Hand Hygiene




### 2. Put on Gown

- Tie neck and waist ties securely




### 5. Put on Gloves

- Put on gloves, taking care not to tear or puncture glove
- If a gown is worn, the glove fits over the gown's cuff



### 3. Put on Mask/N95 Respirator

- Place mask over nose and under chin
- Secure ties, loops or straps
- Mould metal piece to your nose bridge
- For respirators, perform a seal-check



### 4. Put on Protective Eyewear

- Put on eye protection and adjust to fit
- Face shield should fit over brow

